


BIOBASE[®]

Touch Screen Constant-Temperature Incubator



Version 2021

About Us

BIOBASE group, the professional manufacturer for Laboratory and Medical equipment. During long term friendly bilateral cooperation relationship with customers, BIOBASE is always being trusted. The customers purchase BIOBASE made goods, at the mean time, more and more customers are willing to buy a variety of goods from BIOBASE.

In order to satisfy customers demand and ensure customers' benefit, BIOBASE made field trip to many suppliers in China, controlled the product quality strictly in stead of customers, as well as strived from the suppliers to get competitive pricing. Recently, BIOBASE has confirmed many OEM partners and will keep searching for more suppliers upon customers' request.

BIOBASE provide one-stop procurement solution which can not only reduce the resourcing cost, but also can ensure good quality and sufficient after-sale service.

In addition, any abroad strong manufacturer is welcomed to contact us, we can work together to open China market through BIOBASE China.



ISO9001



ISO13485



ISO14001



CE



Touch Screen Constant-Temperature Incubator

Features:

- * Large LCD touch screen.
- * D type: with inner glass door.
G type: with viewing window.
- * Microprocessor controller.
- * Stainless steel chamber.



Double-layer Door (D type)



Glass Viewing Window Door (G type)

Safety Functions:

- * Over temperature alarm.
- * Independent limit temperature alarm system.
- * Leakage current and over-voltage protector.

Standard Accessories:



Test port

Optional Accessories:



Base stand



Universal casters



UV lamp

Technical Parameters:

Model	BJPX-H54BK(D/G)	BJPX-H88BK(D/G)	BJPX-H160BK(D/G)
Capacity	54L	88L	160L
Temperature Control Method	PT100		
Temperature Range	RT+5~85°C		
Temperature Stability	±0.5°C (at 37°C)		
Temperature Uniformity	±1°C (at 37°C)		
Ambient Temperature	5~30°C, recommend 25°C±2°C		
Timing Range	0~999h		
Shelves	2pcs/adjustable		
Power Consumption	345W		495W
Power Supply	AC220V±10%, 50/60Hz(Standard); AC110±10%, 60Hz(Optional)		
Internal Size(W*D*H)mm	395*346*395	500*364*500	558*518*600
External Size(W*D*H)mm	D type: 536*588*681 G type: 522*626*587	D type: 650*573*793 G type: 627*621*692	D type: 718*732*890 G type: 685*779*790
Package Size(W*D*H)mm	D type: 606*648*818 G type: 592*676*724	D type: 720*630*930 G type: 697*671*829	D type: 788*792*1027 G type: 755*829*927
Gross Weight(kg)	D type: 51 G type: 44	D type: 65 G type: 51	D type: 89 G type: 69



C E R T I F I C A T E

ATTESTATION CERTIFICATE OF ELECTROMAGNETIC COMPATIBILITY AND LOW VOLTAGE DIRECTIVES

Technical file of the company mentioned below has been observed
2014/30/EU Electromagnetic Compatibility Directive and
2014/35/EU Low Voltage Directives have been taken as references for these processes

Company Name : **Biobase Bioindustry (Shandong) Co., Ltd.**

Company Address : No.2717 Gongyeyi Road,Mingshui Economic Developing Zone,
Zhangqiu,Jinan,Shandong,China

Related Directives and Annex : **2014/35/EU Low Voltage Directive**
2014/30/EU Electromagnetic Compatibility Directive

Related Standards : **EN 61010-1:2010;EN 61326-1:2013**

Product Name : **Incubator**

Report No and Date : OUCE-SH-174058

Product Brand/Model/Type : BJPX-C50,BJPX-C80,BJPX-C160,BJPX-C270,BJPX-H30,BJPX-H48II,BJPX-H50,
BJPX-H64II,BJPX-H80,BJPX-H123II,BJPX-H160,BJPX-H230II,BJPX-H270,
BJPX-B80,BJPX-B150,BJPX-B200,BJPX-B250,BJPX-B300,BJPX-B400,BJPX-HT150,
BJPX-HT200,BJPX-HT250,BJPX-HT300,BJPX-HT400,BJPX-M150,BJPX-M200,
BJPX-M250,BJPX-M300,BJPX-M400,BJPX-A250,BJPX-A300,BJPX-A400,
BJPX-A500,BJPX-L150,BJPX-L250,BJPX-L300,BJPX-L400,BJPX-L500,BJPX-G,
BJPX-200BX,BJPX-200CX,BJPX-P10,BJPX-P20,BJPX-103B,BJPX-100B,BJPX-200B,
BJPX-1102,BJPX-2102,BJPX-1102C,BJPX-2102C,BJPX-2012,BJPX-2012R

Certificate Number : **M.2016.201.N1827**

Initial Assessment Date : 15.08.2017

Registration Date : 16.08.2017

Reissue Date/No : -

Expiry Date : 15.08.2022


UDEM International Certification
Auditing Training Centre Industry
and Trade Co. Ltd.

The validity of the certificate can be checked through www.udemltd.com.tr. The CE mark shown on the right can only be used under the responsibility of the manufacturer with the completion of EC Declaration of Conformity for all the relevant Directives. This certificate remains the property of UDEM International Certification Auditing Training Centre Industry and Trade Co. Ltd. to whom it must be returned upon request. The above named firm must keep a copy of this certificate for 15 years from the registration of certificate. This certificate only covers the product(s) stated above and UDEM must be noticed in case of any changes on the product(s)

Address: Mutlukent Mahallesi 2073 Sokak (Eski 93 Sokak) No:10 Çankaya - Ankara - TURKEY
Phone: +90 0312 443 03 90 **Fax:** +90 0312 443 03 76
E-mail: info@udemltd.com.tr www.udemltd.com.tr



BIOBASE GROUP

Add: No.9 Gangxing Road, High-tech Zone,
Jinan City, Shandong Province, China
Tel: +86-531-81219801/03 Fax: +86-531-81219804
E-mail: export@biobase.com
Web: www.biobase.cc / www.meihuatrade.com



BIOBASE CHINA



BIOBASE CHINA



BIOBASE CHINA