


**BIOBASE<sup>®</sup>**

## **Disinfection & Sterilization Solution**



Version 2021

# About Us

Founded in 1999, having over 2000 employees (till August, 2018), BIOBASE Group is a new high-tech enterprise focusing on development, production and management of scientific equipment, In Vitro Diagnosis (IVD) instrument and reagent.

BIOBASE Group is specialized in products of 8 areas including medical diagnosis, biosafety protection, disinfection and sterilization solution, water purification system, infant care products, cold chain products, software products, clean room project.

We have full production line for table top steam sterilizer, vertical pressure steam sterilizer, portable pressure steam sterilizer, horizontal pressure steam sterilizer, hot air sterilizer, gas sterilizer, glass bead sterilizer, infrared sterilizer and automatic glassware washer. And every year we have many new products that come into the market.

BIOBASE Group has certificates including ISO9001, ISO13485, ISO14001, SGS, CE, etc for various products.

BIOBASE has set up branch offices in 24 provinces in China, and has also set up warehouses in the United States, India, UAE and Jordan. We now have long-term business partners in more than 130 countries.

BIOBASE Group always adheres to the mission of "service for health for the future" and promotes company's culture vigorously. Providing with high quality products and good service, BIOBASE would like to establish win-win business cooperation with distributors all over the world.



ISO9001



ISO13485



ISO14001



CE



## Hot Air Sterilizer



### Application:

Hot air sterilizer is used to sterilize small surgical instrument, glass, petri dishes etc.

This device is suitable for medical and veterinary clinics, hospitals, pharmacies, health care centers, and laboratories.

### Features:

- \* Inner door without glass, outer door with observation window, back inlet and outlet.
- \* Intelligent microcomputer P.I.D temperature controller, accurate temperature control.
- \* Stainless steel tube circular heater.
- \* The inner box adopts circular arc structure design, and the movable partition can be pulled.

Website: [www.biobase.cc](http://www.biobase.cc) / [www.meihuatrade.com](http://www.meihuatrade.com)

### Technical Parameters:

Model	HAS-T50I	HAS-T70I	HAS-T105I	HAS-T200I
Capacity	48L	80L	105L	180L
Temp. Range	RT+10~250°C			
Temp. Precision	0.1°C			
Temp. Fluctuation	±1°C			
Timing Range	0~9999min			
Shelves No.	2 pcs			
Consumption	1000W	1500W	1600W	2450W
Power Supply	AC220V 50Hz(Standard); 110V 60Hz(Optional)			
Internal Size(W*D*H) mm	419*332*348	450*397*450	550*348*550	602*498*602
External Size(W*D*H) mm	728*595*532	750*680*635	850*620*735	900*775*785
Package Size(W*D*H) mm	810*640*720	840*720*820	940*630*920	990*820*1020
Gross Weight(kg)	60	80	115	145

## **BIOBASE GROUP**

Add: No.9 Gangxing Road, High-tech Zone,  
Jinan City, Shandong Province, China  
Tel: +86-531-81219801/03 Fax: +86-531-81219804  
E-mail: [export@biobase.com](mailto:export@biobase.com)  
Website: [www.biobase.cc](http://www.biobase.cc) / [www.biobase.com](http://www.biobase.com)



BIOBASE CHINA



BIOBASE CHINA



BIOBASE CHINA